QuikRead® CRP Control	Cat. no. 154812	<i>English</i>	154983-1
------------------------------	------------------------	----------------	-----------------

Intended purpose
Intended use: QuikRead® CRP Control is intended for a quality control of quantitative CRP assays by the automated QuikRead go® Instrument or semi-automated QuikRead® 101 Instrument. The control is intended to be used by healthcare professionals in the clinical laboratory and near-patient test (NPT) setting (also known as point-of-care setting). For *in vitro* diagnostic use.

Principles of the procedure: QuikRead CRP Control is with two quantitative instrument-specific assigned values for monitoring the performance of CRP assays used with the QuikRead go Instrument or QuikRead 101 Instrument.
QuikRead go Instrument: Use the control in the same way as an unknown **plasma/serum** specimen would be used, following the instructions for the assay procedure.
QuikRead 101 Instrument: Use the control in the same way as an unknown **whole blood** specimen would be used, following the instructions for the assay procedure.

Reagents

The control is human-based material with added color. It should be handled in the same way as patient sample. The control is stable for 15 minutes in the cuvette.

Content of the kit:

- 1 x 1 ml QuikRead CRP Control

- Instructions for use
The CRP assay and instrument are not provided.

Preparation

Before starting a control assay, read and follow the instructions for use for control, CRP assay and instrument. The control is ready to be used. Allow the material to stabilise for 30 minutes at room temperature (18...25°C). **Mix well.** Do not shake.

Storage

Stable until the expiry date when stored at 2...8°C.

Opened vial: Stable for 2 months at 2...8°C. Do not use any vial that has visible evidence of microbial growth.

QuikRead® CRP Control	Kat. -Nr. 154812	<i>Deutsch</i>
------------------------------	-------------------------	----------------

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Expected values

The individual result for the control should be within the range indi-cated on the control box label. Each laboratory should establish its own mean and precision data. The target values have been set using CRP assays for QuikRead go Instrument and QuikRead 101 Instrument. The values are traceable to ERM®-DA474/IFCC reference material.

Warnings and precautions

- Reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Avoid using the last drops of the control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new vial.
- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- The control contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle the control as if capable of transmitting infectious agents.
- In case of a serious incident, please report it to the manufacturer or its representative and/or national authority.

Disposal

- Disposal of contents according to national and local law.
- Materials of the components
- Paper: Instructions for use
- Cardboard: Outer box
- Glass: Control vials
- Plastic: Vial caps
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

Revision history

Revision history can be found from www.aidian.eu.

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