QuikRead go[®] wrCRP Control High, Cat. no. 153039

English

Intended use

QuikRead go wrCRP Control High is intended for quality control of wrCRP assays by the QuikRead go Instrument. **For *in vitro* diagnostic use.**

Principles of the procedure

Use the control in the same way as unknown **plasma/serum** specimen would be used, following the instructions for the assay procedure.

Reagents

The control is human-based material. It should be handled in the same way as patient sample.

Preparation

QuikRead go wrCRP Control High is ready-to-use. Allow the material to stabilise for 30 minutes at room temperature (18...25°C). **Mix well.** Do not shake.

Storage

Stable until the expiry date when stored at 2...8°C. Opened vial: Stable for 2 months at 2...8°C. Do not use any vial that has visible evidence of microbial growth.

Expected values

The individual result should be within the range indicated on the control vial label. Each laboratory should establish its own mean and precision data.

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153070-2

• Reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.

• Avoid using the last drops of the control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new vial.

• Do not use the product after the expiry date marked on the outer package.

• Do not exceed the stability periods for opened reagents.

• QuikRead go wrCRP Control High contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle QuikRead go wrCRP Control High as if capable of transmitting infectious agents.

Disposal

• Dispose of contents according to national and local law.

• Materials of the components

Paper: Instructions for use

Cardboard: Outer box

Glass: Control vials

Plastic: Vial caps

Rubber: Stopper

• When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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