




QuikRead go[®] CRP




1 Remove the foil cover from the cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette.



2 Fill the blue striped capillary with the sample (20 µL) up to the white stopper. Ensure that there are no air bubbles in the capillary. Wipe away any residual sample from the outside of the capillary.




3 Place the sample (20 µL) into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.




4 Close the cuvette tightly with the CRP reagent cap. Do not press down the blue inner part of the reagent cap. Perform the test within 2 hours.



5 Choose **Measure** on the display of the QuikRead go instrument.



6 Insert the cuvette into the measurement well of the instrument. The barcode should be facing you.



7 The result appears on the display and the cuvette rises up automatically when the measurement is completed.

Please consult the instructions for use before performing the test.

QuikRead go[®] CRP



Kit

Unopened kit at 36–46°F (2–8°C): until the **expiry date** of the kit



Reagent cap tube

Storage at 36–46°F (2–8°C): until the **expiry date** of the kit
Storage at 64–77°F (18–25°C): **1 month** (3 months when 7.5 hours/day)



Cuvettes

Unopened foil pouch at 36–77°F (2–25°C): until the **expiry date** of the kit
Unopened cuvettes after opening the foil pouch: at 36–46°F (2–8°C): **6 months** / at 64–77°F (18–25°C): **3 months**



Opened cuvette

Let the buffer warm up to 64–77°F (18–25°C) before testing
Perform the test within **2 hours** from opening the cuvette
The sample is stable in the buffer for **2 hours**



QuikRead go CRP Control Set

Unopened vial at 36–46°F (2–8°C): until the **expiry date**
Opened vial at 36–46°F (2–8°C): up to **1 month**



QuikRead go CRP Verification Set

Unopened vial at 36–46°F (2–8°C): until the **expiry date**
Opened vial at 36–46°F (2–8°C): up to **3 weeks**

Contact info

Distributed in USA by
Clinical Diagnostic Solutions, Inc.

1800 NW 65th Avenue, Plantation, FL 33313

For information contact: sales@boule.com
www.csdolinc.com

Product	Cat. no.
QuikRead go CRP with capillaries and plungers	145215
QuikRead go CRP Control Set includes 2 CRP control levels	145217
QuikRead go CRP Verification Set includes 3 CRP levels for calibration verification	145216
QuikRead go Instrument	145218