Hygicult[®] Y&F (Easicult M)

Instructions for use • English

Hygicult[®] Y&F (Easicult M)

Intended use

Hygicult Y&F (Easicult M) slides are intended for rapid monitoring of fungal contamination in different types of materials, both solid and liquid. The test can be performed on-site, or the slides can be used as convenient transport media for samples.

The slide is covered on both sides with Malt Agar which supports rapid growth of yeasts and moulds. The growth of bacteria is inhibited. The main significance of the test is that elevation of levels of fungi can be detected. Normal levels must therefore be first established.

Contents of the kit

Hygicult Y&F (Easicult M)	Cat. No. 68013
Test slides	10 pcs
Labels	10 pcs
Instructions for use	1 pc

Typical formulation

Malt agar	
Malt agar	Antibiotics
Yeast extract	Agar agar
Glucose D	Water
Lactic acid	

Warnings and precautions

Do not use product beyond the expiry date marked on the kit. Do not use the kit if you notice

- discoloration or dehydration of the growth medium
 - detachment of the growth media from the plastic slide
 - evidence of bacterial or fungal growth
- Do not touch the growth because any colony growing on the slide may be pathogenic.

Storage

Store the kit at room temperature (18...25°C / 64...77°F) protected from draught, temperature fluctuations and light sources.Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box and on the cap of each slide.

Sampling

To avoid contamination, the growth medium should not come into contact with any other material than the one to be tested. On the other hand, it is important that the growth medium makes full contact with the material to be tested. After sampling screw the slide tightly back into the tube.

Contact inoculation (Fig. 1a, 1b)

Solid surfaces can be tested by pressing each side of the slide firmly against the surface for three or four seconds. The slide should be held still during pressing. The hinged design offers ease of use.

Dipping (Fig. 2)

Fluid samples are tested by dipping the slide in the liquid for three or four seconds. Blot the last drops on absorbent paper.

Swabbing (Fig. 3)

Semisolid materials or objects that are difficult to reach can be tested by carefully rolling a sterile swab over an area delimited using e.g. a frame. If the object is dry, the swab should first be moistened with sterile water. The moistened swab can also be used for obtaining samples from powders (e.g. spices) or viscous fluids.

After swabbing the sample area, roll the swab gently over the agar surfaces of the slide from left to right and from bottom to top.

Incubation (Fig. 4)

Incubate the slide tightly enclosed in its tube at 27...30°C. Higher temperatures are not recommended since some of the most common yeasts and moulds do not grow at high temperatures. The incubation time is three to five days. Growth is usually visible within three days but a five-day incubation is recommended to detect any slow-growing fungi.

Interpretation of results (Fig. 5)

Remove the slide from its tube after incubation. The slide must not be shaken since any spores produced by the fungi can easily spread yielding falsely high fungal counts. Determine the fungal count (number of colony forming units, CFU) by comparing the density of growth on the slide with the model chart. No bacterial colonies should be visible, as the bacterial growth is inhibited by low pH and the presence of antibiotics on the slide. The growth appearing on the Hygicult Y&F (Easicult M) slide may consist of mould only, yeast only or a mixture of the two. Mould colonies are soft and fluffy and usually pale, green or black in colour. Yeasts usually grow in dome-shaped colonies but may sometimes be flat and dry. The colonies are often pale or red.

Limitations of the method

Results obtained with different inoculation systems should not be compared. Valid comparisons can only be made among results obtained using the same technique on the same type of material.

Disposal

- Dispose of contents according to national and local law.
- All used components should be handled and disposed of as potentially pathogenic material.
- Materials of the components:

Paper: Instructions for use, patient labels

Cardboard: Kit box

Plastic: Tubes, caps and diplslides
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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Hygicult® is a registered trademark of Aidian Oy.
Aidian Oy