

AIDIAN



Cultura M

A small and handy incubator for dipslides and contact slides.

A small and handy incubator for dipslides and contact slides.

Cultura M, a small, durable and reliable incubator is suitable for the incubation of Aidian's contact slide Hygicult® and dip slides Uricult®, and Easicult®.

The transparent door allows you to check inside the incubator without opening the door. Cultura M is small and light and can be easily transported wherever you need to use it.

Cultura M is equipped with a thermometer and multirack shelf, which are available also separately.



Cultura M incubator

Temperature range	25...45°C
Temperature stability	± 1°C
Voltage	220-240 V / 100-120 V
Frequency	50/60 Hz
Voltage	30W
Safety class	II
Dimensions (mm)	external W: 310 H: 155 D: 168
	internal W: 220 H: 120 D: 150
Case	Flame resistant PBDE-free PC & ABS
Door	PC, transparent
Thermometer	up to 60°C, in special tube
Base cover	ABS
Warranty	2 years from initial installation
Weight	1.1 kg



This CE Mark indicates that this equipment complies with the following directive:

- EMC 2014/30/EU
- LVD 2014/35/EU
- RoHS III 2015/863/EU
- WEEE II 2012/19/EU
- EU REACH 1907/2006
- UK SI 2016 No. 1091
- UK SI 2016 No. 1101

Contact info

Aidian Oy

aidian@aidian.eu / +358 10 309 3000

For orders: orders@aidian.eu

Product	Cat. no	Order no
Cultura M, 220-240 V	05866	129997
Cultura M, 100-120 V	06266	132754
Cultura M, 220-240 V with UK plug	-	133192
Thermometer, accessory	05724	129981
Multirack, accessory	-	139248

AIDIAN

aidian.eu