


Get ready for flu season!



AIDIAN

Tools for the correct diagnosis.

As the symptoms and signs of influenza, COVID-19 and other respiratory infections somewhat overlap, the correct diagnosis of the infection is crucial.

Aidian offers a range of respiratory panels for the molecular diagnostics of influenza, COVID-19 and other respiratory infections. ELITech's ELITe InGenius® and BeGenius® platforms together with SARS-CoV-2 PLUS (SARS-CoV-2, Flu A, Flu B, RSV) and the Respiratory Viral PLUS (Flu A, Flu B, RSV, hMPV) ELITE MGB® panels offer a fully automated sample to result solution. Mikrogen's *ampliCube* Respiratory Flu & SARS-CoV-2 assay is a flexible solution for differentiation of SARS-CoV-2, Flu A and Flu B with one test. Mikrogen's Pathogen group Respiratory F contains a set of panels for the detection of SARS-CoV-2, Flu A, Flu B and RSV. All ELITech's and Mikrogen's respiratory panels are fully compatible with open platforms.

SD BIOSENSOR's STANDARD™ M10 is an automated point of care random access system, that has both real-time PCR and isothermal amplification in one platform. STANDARD M10 FLU/RSV/SARS-COV-2 is a fast multiplex PCR assay for differentiation of SARS-CoV-2, RSV, Flu A and Flu B in one test. A swift diagnosis ensures the right treatment and minimal spread of the infection.

Respiratory panels	Product no.	Order. no.
ELITech ELITe InGenius	INT030	153717
ELITech ELITe BeGenius	INT040	155183
ELITech SARS-CoV-2 PLUS ELITe MGB Kit	RTS180ING	155028
ELITech Respiratory Viral PLUS ELITE MGB Kit	RTS160ING	155030
Mikrogen <i>ampliCube</i> Respiratory Flu & SARS-COV-2	50163 / 50173 (LC)	155033 / 155034
Mikrogen Pathogen group Respiratory F (including RSV)	59976 (50143 & 50152)	155035
SD BIOSENSOR STANDARD M10 Console	11M1011	155067
SD BIOSENSOR STANDARD M10 Module	11M1012	155068
SD BIOSENSOR M10 FLU/RSV/SARS-COV-2	11FLU10A	155176
SARS-CoV-2 tests	Product no.	Order. no.
ELITech SARS-CoV-2 ELITe MGB Kit	RTS170ING-96	154715
COVID-19 Plus Real Amp Kit	IFMR-45	154600
Mikrogen <i>ampliCube</i> Coronavirus SARS-CoV-2	50143	154612
SD BIOSENSOR M10 SARS-COV-2	11COV10A	155069

Aidian Oy (Finland)
info@aidian.fi / 050 380 9684
tilaukset@aidian.fi

Aidian Sweden AB
info@aidian.se / 08 623 0400
order@aidian.se

Aidian Norway AS
info@aidian.no / 066 78 5630
ordre@aidian.no

Aidian Denmark ApS
info@aidian.dk / 08 610 3020
order@aidian.dk