



# Get ready for flu season!



**AIDIAN**

# Tools for the correct diagnosis.

As the symptoms and signs of influenza, COVID-19 and other respiratory infections somewhat overlap, the correct diagnosis of the infection is critical.

Aidian offers a range of respiratory panels for the molecular diagnostics of influenza, COVID-19 and other respiratory infections. ELITech's ELITe InGenius® and BeGenius® platforms together with SARS-CoV-2 PLUS (SARS-CoV-2, Flu A, Flu B, RSV) and the Respiratory Viral PLUS (Flu A, Flu B, RSV, hMPV) Elite MGB® panels offer a fully automated sample to result solution. Mikrogen's *ampliCube* Respiratory Flu & SARS-CoV-2 assay is a flexible solution for differentiation of SARS-CoV-2, Flu A and Flu B with one test. Mikrogen's Pathogen group Respiratory F contains a set of panels for the detection of SARS-CoV-2, Flu A, Flu B and RSV. All ELITech's and Mikrogen's respiratory panels are fully compatible with open platforms.

SD BIOSENSOR's STANDARD™ M10 is an automated point of care random access system, that has both real-time PCR and isothermal amplification in one platform. STANDARD M10 FLU/RSV/SARS-COV-2 is a fast multiplex PCR assay for differentiation of SARS-CoV-2, RSV, Flu A and Flu B in one test. Qlife's egoo.health system uses Aidian's proprietary SIBA® technology and disposable Egoo capsules for near patient molecular analysis of SARS-CoV-2 in 30 minutes with less than two minutes hands-on time.

A swift diagnosis ensures the right treatment and minimal spread of the infection.

Respiratory panels	Product no.	Order. no.
ELITech ELITe InGenius	INT030	153717
ELITech ELITe BeGenius	INT040	155183
ELITech SARS-CoV-2 PLUS ELITe MGB Kit	RTS180ING	155028
ELITech Respiratory Viral PLUS ELITe MGB Kit	RTS160ING	155030
Mikrogen <i>ampliCUBE</i> Respiratory Flu & SARS-COV-2	50163 / 50173 (LC)	155033 / 155034
Mikrogen Pathogen group Respiratory F (including RSV)	59976 (50143 & 50152)	155035
SD BIOSENSOR STANDARD M10 Console	11M1011	155067
SD BIOSENSOR STANDARD M10 Module	11M1012	155068
SD BIOSENSOR M10 FLU/RSV/SARS-COV-2,10 tests	11FLU10A	155176
SARS-CoV-2 tests	Product no.	Order. no.
ELITech SARS-CoV-2 ELITe MGB Kit	RTS170ING-96	154715
COVID-19 Plus Real Amp Kit	IFMR-45	154600
ELITech SARS-CoV-2 Variants ELITe MGB Kit	RTS171ING	154833
ELITech SARS-CoV-2 Extended ELITe MGB Kit	RTS172ING	154987
Mikrogen <i>ampliCube</i> Coronavirus SARS-CoV-2	50143	154612
Mikrogen <i>ampliCube</i> SARS-COV-2 Variants	50145	155039
Mikrogen <i>ampliMelt</i> SARS-COV-2 Variants 1	50147	155040
Mikrogen <i>ampliMelt</i> SARS-COV-2 Variants 2	50141	155041
SD BIOSENSOR M10 SARS-COV-2, 10 tests	11COV10A	155069
Qlife Egoo instrument	E78852	154845
Qlife SARS-CoV-2 capsule, 10 tests	E48364	154843
Qlife SARS-CoV-2 Lysis buffer, 50 pc	E70464	154844

Aidian Oy (Finland): [info@aidian.fi](mailto:info@aidian.fi) / 050 380 9684 / For orders: [tilaukset@aidian.fi](mailto:tilaukset@aidian.fi)

Aidian Sweden AB: [info@aidian.se](mailto:info@aidian.se) / 08 623 04 00 / For orders: [order@aidian.se](mailto:order@aidian.se)

Aidian Norway AS: [info@aidian.no](mailto:info@aidian.no) / 066 78 56 30 / For orders: [ordre@aidian.no](mailto:ordre@aidian.no)

**AIDIAN**