

AIDIAN



QuikRead go[®] CRP

C-Reactive Protein, a Proven Tool in Infection Diagnosis

An Easy to Use Test for Quantitative Measurement of CRP.


C-reactive protein (CRP) is a well-known and used biomarker for evaluating inflammation, infection and tissue injury¹. It is a fast and reliable reflector of the acute-phase response to the treatment². CRP is normally present at very low concentrations in the blood of healthy people³. Concentrations are markedly increased in acute inflammatory processes⁴. CRP supports in reaching the correct diagnosis and treatment decision^{1,5}. CRP provides information on the course of bacterial infections and is a good indicator of the efficacy of antimicrobial therapy^{2,6}.



QuikRead go CRP is Reliable, Fast and Easy to Use

- Proven quality – substantially equivalent to two clinical chemistry analyzers⁷
- Reproducible quantitative CRP result within range of 5–150 mg/L⁷
- Performed on a small blood sample (whole blood, plasma or serum)⁷
- Test time 2 minutes
- Minimum hands-on time
- Easy testing procedure, intuitive user interface and maintenance-free instrument⁷
- Bi-directional connectivity to LIS/HIS⁷

For more detailed information, please consult the instructions for use.



Sample



Measure



Result

References

1. Aabenhus R et al. Cochrane Database Syst Rev. 2014;11.
2. van Leeuwen MA et al. Baillieres Clin Rheumatol. 1994;8(3):531–52.
3. Shine B et al. Clin Chem Acta. 1981;117(1):13–23.
4. Pebys MB. Oxford Textbook of Medicine. 2003; p. 150–156.
5. Little P et al. Lancet. 2013;382(9899):1175–82.
6. Olaison L et al. Arch Intern Med. 1997;157(8):885–92.
7. QuikRead go CRP instructions for use.

Contact info

Aidian Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83
FI-02101 Espoo, Finland
aidian@aidian.eu
+358 10 309 3000
www.aidian.eu

Distributed in USA by

Clinical Diagnostic Solutions, Inc.
1800 NW 65th Ave, Plantation
FL 33313
+1 954-791-1773
sales@boule.com
www.cdsolinc.com

Product	Cat. no
QuikRead go CRP, 50 tests	145215
QuikRead go CRP Control Set	145217
QuikRead go CRP Verification Set	145216
QuikRead go Instrument	145218

QuikRead go CRP test is cleared to be used by trained professionals in clinical laboratories. CLIA categorization: moderate complexity. CRP results must be evaluated in the light of clinical information of the patient.

AIDIAN