


Certificate

Aidian Oy

using

QuikRead go HbA1c

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.



Total Error =	3.3	mmol/mol
Bias =	-0.8	mmol/mol
Imprecision =	2.5	%
Grade =	Bronze	

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61 : 752-59)

Date of Certification : 01 June 2022

Date of Expiry : 01 June 2023


IFCC Network Coordinator
C. Siebelder