


## 2019-nCoV IgG/IgM Rapid Test Cassette

CE-IVD  
Marked

To support in rapid diagnosis of the new coronavirus pneumonia.

The 2019-nCoV IgG/IgM rapid test offers high accuracy and convenient usability as only a fingerstick blood sample is needed.

Results in only 10 minutes.




### Product information

Sample: Fingerstick blood 20µl

Accuracy compared to PCR

- IgG sensitivity 100%, specificity 98%
- IgM sensitivity 85%, specificity 96%

Reference: Instructions for use.



### General information about IgG and IgM in the coronavirus infection

Antibody detection was done for 16 patients that were positive for 2019-nCoV in molecular test:

- IgM and IgG titers were relatively low or undetectable on the first day of sampling. By day 5, the titers of both antibodies had risen in most patients.
- The rate of detection rose from 50% (8/16) on day 0 to 81% (13/16) on day 5 for IgM, and from 81% (13/16) to 100% (16/16) for IgG.

Reference: Medscape: Novel Coronavirus: Study Suggests Multiple Shedding Routes  
Nicola M. Parry, DVM Feb 21,2020 <https://doi.org/10.1080/22221751.2020.1729071>

Please contact your local Aidian office for further inquiries

Finland: [info@aidian.fi](mailto:info@aidian.fi) tel. 050 380 9684 / Sweden: [info@aidian.se](mailto:info@aidian.se) tel. 08 623 04 00  
Norway: [info@aidian.no](mailto:info@aidian.no) tel. 66 78 56 30 / Denmark: [info@aidian.dk](mailto:info@aidian.dk) tel. 86 10 30 20

# AIDIAN

## Real time PCR assay for detection of the 2019-nCoV

CE-IVD  
Marked


Compatible with all the main Real Time PCR platforms

Sample material for RNA extraction: Respiratory specimens  
such as alveolar lavage fluid, throat swab, sputum

The PCR kit contains internal control, positive and negative controls

Processable on the ELITE InGenius - PCR automation

- Sample to result
- 12 samples in one run in 2.5 hours  
with only minimal hands on steps
- Universal extraction for further analysis  
of the RNA/DNA eluate



Please contact your local Aidian office for further inquiries

Finland: [info@aidian.fi](mailto:info@aidian.fi) tel. 050 380 9684 / Sweden: [info@aidian.se](mailto:info@aidian.se) tel. 08 623 04 00  
Norway: [info@aidian.no](mailto:info@aidian.no) tel. 66 78 56 30 / Denmark: [info@aidian.dk](mailto:info@aidian.dk) tel. 86 10 30 20