

QuikRead go  
easy CRP



# QuikRead go<sup>®</sup> easy CRP

Test before prescribing

**ORION**  
DIAGNOSTICA

# Test before prescribing

QuikRead go easy CRP is an easy to use and reliable point-of-care test, which gives quantitative CRP results in 2 minutes. CRP provides information for the detection and evaluation of infection and inflammatory disorders. Measurement of CRP can also be used to promote more prudent use of antibiotics and to support physician to reach the right diagnosis and correct treatment decision during patient consultation.

## Because...

### UNIQUE ONE-STEP SAMPLE COLLECTOR

- » Fast and easy one-step sample collection improves workflow
- » Convenient for the patient
- » Flexible testing as the sample can be kept in the cuvette for two hours before measurement

### ROOM TEMPERATURE STORAGE

- » Easy to store in the health care unit
- » Enables immediate use, no need to warm up the reagents

### RELIABLE CRP RESULT IN 2 MINUTES

- » CRP measuring range from whole blood is 1–200 mg/l and 1–120 mg/l with serum or plasma samples

### QUIKREAD GO MULTIANALYTE POINT-OF-CARE SYSTEM

- » Portable and fully automatic point-of-care system
- » Bi-directional connectivity to most HIS and LIS systems



**QuikRead go**<sup>®</sup>

## EASY ONE-STEP SAMPLING



Sample



Measurement



Result

*For more detailed information, please consult the instructions for use.*

Product	Cat. no.
QuikRead go easy CRP, 50 tests	
» With 10 µl Sample Collector	153287
Ultrasensitive CRP Control	68257
QuikRead go wrCRP Control	153038
QuikRead go wrCRP Control High	153039
QuikRead go CRP Control Low	153765
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763
QuikRead go Instrument	133893



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

QuikRead go<sup>®</sup> is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy  
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland  
[www.oriondiagnostica.com](http://www.oriondiagnostica.com), [www.quikread.com](http://www.quikread.com), [www.tacklearm.com](http://www.tacklearm.com)